






YPB-2100

VISUAL CHART

YEOSN



Color Blindness Test

YPB-2100 offers more than 70 various charts for color blindness test. All color blindness charts are accompanied with easy-to-use instructions.

various Types of test Charts

YPB-2100 presents 150 various types of charts, including Snellen E, Landolt C, Letter, Number, children chart, standard ETDRS and multi-function charts. High resolution LCD enables wide range of tests.

Specification

Display	23.8", 1920×1080
Max Brightness	300cd/m ²
Working distance	1.5m ~ 6m, step 0.1m
Fixation	Wall-mounted, Desktop(Optional)
Standard Chart	Snellen E, Landolt C, Letter, Numbers, Children chart, ETDRS
Functional Chart	Maddox Rod, Cross Cylinder, Worth 4, Stereoscopic, etc.
Display type	Arrow, Single, Horizontal, Vertical, Full, Random
Chart Mask	R/G, B/W, Low contrast, Low illumination
Dimension (Host)	595.6(L) 388.7(W) 62(H) mm/4.65kg
Dimension (Remote controller)	186(L) 55(W) 17.2(H) mm/80g
Power supply	AC ~ 100-240V, 50/60Hz, Power 1.0-0.5A



U2, F3, No.33, 30th Tir St, Tehran
www.binateb.com

Shop@binateb.com
(021)-68 25 8000

