

LIANYUNGANG

TIANNNUO

OPTICAL INSTRUMENT CO., LTD.



TianNuo Optical Instrument



Trial lens Set PL-232-BI

* plastic mounts BI with power marks branded clearly on both sides of each handle



Trial lens Set PL-232-BV

* Nylon mounts BV with power marks branded clearly on both sides of each handle, Cylinder and Prismatic lens in square shape to prevent rotation

Specification

* Reduced-aperture trial lenses of 22mm in diameter, 232 pieces included, of accurate power degrees and perfect outlook finish

* Standard removable tray with grey velvet and labels of silver background

* Standard carrying case ALCC No,009 of aluminum framework, with black ABS cover, black inner-felt and combination locks

* Dimensions: Tray: W 518mm* L 318mm* H 35mm Case: W 545mm* L 345mm* H 95mm * Weight: 5.1kg

Concave(-) Sphere				Convex(+) Sphere				concave(-) Cylinder		Convex(+) Cylinder		Prism		Accessory	
List	Pcs	List	Pcs	List	Pcs	List	Pcs	List	Pcs	List	Pcs	List	Pcs	Specification	Pcs
0.12	2	6.00	2	0.12	2	6.00	2	0.12	2	0.12	2	1.00	1	Occluder	2
0.25	2	6.50	2	0.25	2	6.50	2	0.25	2	0.25	2	2.00	1	Pin Hole (PH1.5mm)	1
0.50	2	7.00	2	0.50	2	7.00	2	0.50	2	0.50	2	3.00	1	Pin Hole (PH1.5mm)	1
0.75	2	8.00	2	0.75	2	6.00	2	0.75	2	0.75	2	4.00	1	Slit (SS1.0 Or SS1.5)	1
1.00	2	9.00	2	1.00	2	9.00	2	1.00	2	1.00	2	5.00	1	Red Filter Lens	1
1.25	2	10.00	2	1.25	2	10.00	2	1.25	2	1.25	2	6.00	1	Green Filter Lens	1
1.50	2	11.00	2	1.50	2	11.00	2	1.50	2	1.50	2	7.00	1	White Maddox	1
1.75	2	12.00	2	1.75	2	12.00	2	1.75	2	1.75	2	8.00	1	plane lens	1
2.00	2	13.00	2	2.00	2	13.00	2	2.00	2	2.00	2	9.00	1	Frosted Lens	1
2.25	2	14.00	2	2.25	2	14.00	2	2.25	2	2.25	2	10.00			
2.50	2	15.00	2	2.50	2	15.00	2	2.50	2	2.50	2				
2.75	2	16.00	2	2.75	2	16.00	2	2.75	2	2.75	2				
3.00	2	18.00	2	3.00	2	18.00	2	3.00	2	3.00	2				
3.25	2	20.00	2	3.25	2	20.00	2	3.25	2	3.25	2				
3.50	2			3.50	2			3.50	2	3.50	2				
3.75	2			3.75	2			4.00	2	4.00	2				
4.00	2			4.00	2			4.50	2	4.50	2				
4.50	2			4.50	2			5.00	2	5.00	2				
5.00	2			5.00	2			6.00	2	6.00	2				
5.50	2			5.50	2										
Total: 68 Pcs				Total: 68 Pcs				Total: 38 Pcs		Total: 38 Pcs		Total: 10 Pcs		Total: 10 Pcs	



U2, F3, No.33, 30th Tir St, Tehran

www.binateb.com

Shop@binateb.com

(021)-68 25 8000





Trial Clip TC1

- * With proper pinchcock on lens holders
- * A twin of clips supplied simultaneously
- * Available for 3 pieces trial lenses to be inserted at each side



Prism Bar HB16&VB15


In a rigid protective case, complete with:

- 1 vertical bar with the following prism diopters: 1 Δ , 2 Δ , 3 Δ , 4 Δ , 5 Δ , 6 Δ , 8 Δ , 10 Δ , 12 Δ , 14 Δ , 16 Δ , 18 Δ , 20 Δ , 25 Δ , 30 Δ .
- 1 horizontal bar with the following prism diopters: 1 Δ , 2 Δ , 4 Δ , 6 Δ , 8 Δ , 10 Δ , 12 Δ , 14 Δ , 16 Δ , 18 Δ , 20 Δ , 25 Δ , 30 Δ , 35 Δ , 40 Δ , 45 Δ .



UTF5080 Pro

- * PD scale: 50mm ~ 80mm
- * Net weight : 55g
- * knobs fastened with screws
- * Back lens holders of round corners
- * Easily-shaped ear rests
- * Available for 5 pieces trial lenses to be inserted at each side



UTF5080

- * PD scale: 50mm ~ 80mm
- * Net weight : 61g
- * knobs fastened with screws
- * Back lens holders of round corners
- * Easily-shaped ear rests
- * Available for 5 pieces trial lenses to be inserted at each side

